

NordVal Certificate

Issued for:	TRANSIA TM PLATE <i>Salmonella</i> Gold
NordVal No:	001
First approval date:	14 June 2001
Renewal date:	1 June 2007
Valid until:	1 June 2009

Oslo, 1. June 2007



Raisio Diagnostics AB, FO Petersons Gata 32m SE 421 31 Västra Frölunda, Sweden, has applied for renewal of the NordVal certificate for the TRANSIATM PLATE *Salmonella* Gold. The TRANSIATM PLATE *Salmonella* Gold is produced by Raisio Diagnostics SAS, 8 rue Saint Jean de Dieu, 69007 Lyon, France.

NordVal has studied the application for renewal as well as the enclosures. NordVal has concluded that it is demonstrated that the requirements of the NordVal validation protocol are fulfilled. The TRANSIATM PLATE *Salmonella* Gold was tested against ISO 6579:2002: Horizontal method for the detection of *Salmonella* spp.

The statistical interpretation of the comparisons of these methods was conducted by AFNOR in 2001 and 2005, respectively. No statistical difference in the results between ISO 6579:2002 and the TRANSIATM PLATE *Salmonella* Gold was found. The limit of detection: 0.3 -1.4 cfu/25g, the specificity rate about 100% and the sensitivity rate >95%. The method is applicable to food, feed and environmental samples.

Please note that NordVal should be informed of any changes of the method. Further that the certificate is valid for two years.

Yours sincerely


Sven Qvist
Chair of NordVal
Hilde Skår Norli
NMKL Secretary General